

# MycoClean™ Mycoplasma Prevention Spray

MycoClean Mycoplasma Prevention Spray is the lab's solution for decontamination.

It can be used for cell culture areas, incubators, work benches, storage containers and apparatuses.



- Ready-to-use / Convenient and Stable
- Non-corrosive and Non-carcinogenic
- Can be stored at RT for 2 year
- Manufactured using safe ingredients
- Elimination range: Mycoplasma, bacteria, fungi, viruses

**MycoClean™ Mycoplasma Prevention Spray** is an effective solution for reliable disinfection of laboratory surfaces and apparatus including clean benches, incubators, work benches, cell storage boxes and liquid nitrogen containers. MycoClean™ shows an excellent effectiveness against mycoplasma. Contamination-free and clean surfaces can be achieved in a single application.

- Mycoplasma lacks a cell wall, making them unresponsive to most antibiotic treatments
- Prevention spray forms pores in cell membrane using osmotic pressure and eradicate mycoplasma
- Using with this easy and convenient spray product, contamination-free surfaces can be achieved with a single application

A test using the **6 mycoplasma species** after multiplying them on the plate and compared with 70% EtOH (For each dish, The right half of the area was sprayed out by MycoClean). As a result, the area **sprayed by MycoClean does not show any colony of mycoplasma**. This test also have demonstrated that 70% EtOH which is generally used for removing is not effective to mycoplasma.

Product	Cat. #	Qty.
MycoClean™ Mycoplasma Prevention Spray	21083	200mL × 1 ea
e-Myco™ Mycoplasma PCR Detection Kit	25235	48tubes
e-Myco™ plus Mycoplasma PCR Detection Kit	25237	48tubes
e-Myco™ VALiD Mycoplasma PCR Detection Kit	25239	48tests
e-Myco™ VALiD-Q Mycoplasma qPCR Detection Kit	25245	50test

## How to order

iNtRON Biotechnology, Inc | T. +82 31.739.5077 | Fax : +83 31.739.5264 | [tech@intronbio.com](mailto:tech@intronbio.com)